



[www.ESECO-Speedmaster.com](http://www.ESECO-Speedmaster.com)

ONE ESECO ROAD  
CUSHING, OKLAHOMA  
74023-9912 U.S.A.  
PHONE: 918-225-1266  
FAX: 918-225-1284

E-MAIL: [info@eseco-speedmaster.com](mailto:info@eseco-speedmaster.com)

## ST-21

### 21-STEP TRANSMISSION DENSITY TABLET

Cat. No. 5004 \$145.00



1	0.15
2	0.39
3	0.60
4	0.73
5	0.91
6	1.12
7	1.30
8	1.47
9	1.61
10	1.83
11	2.00
12	2.24
13	2.42
14	2.62
15	2.84
16	3.05
17	3.33
18	3.51
19	3.72
20	4.12
21	4.72

011216-1

**CALIBRATION TRACEABLE TO  
U.S. NATIONAL INSTITUTE OF STANDARDS AND TECHNOLOGY**

The Density Tablet is an 8" X 1" X-Ray Film Strip divided into 21 progressive density steps. The value of the least dense step is in the range of .15 and the value of the most dense step is in the range of 5.00.

One of the more usual applications of the ST-21 is its use as a Certified Standard Reference Material for verification of the accuracy of Densitometers and associated Density Calibration Check Tablets.